

ANORGANIKUS FOSZFOR STABIL FOLYÉKONY  
STANDARD  
INORGANIC PHOSPHORUS STABLE LIQUID  
STANDARD



Kódszám/Catalog No.: 52101  
Kiszerezés/Size: 1x5 ml

Gyártási szám: 7135 #####

Lejáratidő: ÉÉÉÉ-HH

**Felhasználás:**

Anorganikus foszfor reagens kalibrálására szolgáló stabil folyékony foszfor standard.  
Felhasználásra kész.  
A standard koncentrációja LOT specifikus.

**Koncentráció:**

Reagens	Érték	Egység
Anorganikus foszfor stabil folyékony reagens (UV)	###	mmol/l

**Visszavezethetőség:**

Belső standard


**Stabilitás:**


Felbontás nélkül 2-8°C-on tárolva a lejáratidőig.  
Felbontás után 2-8°C-on tárolva minimum 4 hétig.


**Figyelmeztetések:**

Csak in vitro diagnosztikai célokra használható.  
Ne használja a standardet a lejáratidő eltelte után!

**A címkén a következő szimbólumok találhatóak:**

 In vitro diagnosztikai használatra

 Lejáratidő (a hónap utolsó napja)

 Tárolási hőmérséklet

 Gyártási szám

 Kódszám

Lot No: 7135 #####

Expiry date: YYYY-MM

**Preparation:**

Phosphorus standard for the calibration of inorganic phosphorus reagent.  
Ready to use.  
Concentration is LOT specific.

**Concentration:**

Reagent	Value	Unit
Inorganic phosphorus stable liquid reagent (Molybdate UV)	###	mg/dl

**Traceability:**

Internal standard


**Stability:**


Without opening at 2-8°C until the expiry date.  
After first opening at least 4 weeks at 2-8°C.


**Warnings:**

For in vitro diagnostic use only.  
Do not use after the expiry date!

**The following symbols are used on labels:**

 For in vitro diagnostic use

 Expiry date (last day of the month)

 Temperature limitation

 Batch Code

 Code