

Gyártási szám: 7685 #####

Lejáratidő: ÉÉÉÉ-HH

**Felhasználás:**

Glükóz reagens kalibrálására szolgáló stabil folyékony glükóz standard.

Felhasználásra kész.

A standard koncentrációja LOT specifikus.

**Koncentráció:**

Reagens	Érték	Egység
Glükóz stabil folyékony reagens (GOD/PAP)	##	mmol/l
Glükóz stabil folyékony reagens (HK/G6PDH)	##	mmol/l

**Visszavezethetőség:**

IDMS

**Stabilitás:**


Felbontás nélkül 2-8°C-on tárolva a lejáratidőig.


Felbontás után 2-8°C-on tárolva minimum 4 hétig.


**Figyelmeztetések:**

Csak in vitro diagnosztikai célokra használható.


Ne használja a standardet a lejáratidő eltelte után!

**A címkén a következő szimbólumok találhatóak:**
 In vitro diagnosztikai használatra

 Lejáratidő (a hónap utolsó napja)

 Tárolási hőmérséklet

 Gyártási szám

 Kódszám

Lot No: 7685 #####

Expiry date: YYYY-MM

**Preparation:**

Glucose standard for the calibration of glucose reagent.

Ready to use.

Concentration is LOT specific.

**Concentration:**

Reagent	Value	Unit
Glucose stable liquid reagent (GOD/PAP)	#	mg/dl
Glucose stable liquid reagent (HK/G6PDH)	#	mg/dl

**Traceability:**

IDMS

**Stability:**


Without opening at 2-8°C until the expiry date.


After first opening at least 4 weeks at 2-8°C.

**Warnings:**

For in vitro diagnostic use only.

Do not use after the expiry date!

**The following symbols are used on labels:**
 For in vitro diagnostic use

 Expiry date (last day of the month)

 Temperature limitation

 Batch Code

 Code