

Kódszám/Catalog No.: 50511
Kiszerezés/Size: 1x5 ml

Gyártási szám: 7675 #####

Lejárati idő: ÉÉÉÉ-HH

Felhasználás:

Húgysav reagens kalibrálására szolgáló stabil folyékony húgysav standard.

Felhasználásra kész.

A standard koncentrációja LOT specifikus.

Koncentráció:

Reagens	Érték	Egység
Húgysav stabil folyékony reagens (ADPS)	#	μmol/l

Visszavezethetőség:

GC-IDMS

Stabilitás:


Felbontás nélkül 2-8°C-on tárolva a lejárató időig.
Felbontás után 2-8°C-on tárolva minimum 4 hétig.


Figyelmeztetések:


Csak in vitro diagnosztikai célokra használható.

Ne használja a standardet a lejárató dátum eltelte után!

A címkén a következő szimbólumok találhatóak:

 In vitro diagnosztikai használatra

 Lejárató idő (a hónap utolsó napja)

 Tárolási hőmérséklet

 Gyártási szám

 Kódszám

Lot No: 7675 #####

Expiry date: YYYY-MM

Preparation:

Uric acid standard for the calibration of uric acid reagent.

Ready to use.

Concentration is LOT specific.

Concentration:

Reagent	Value	Unit
Uric acid stable liquid reagent (ADPS)	###	mg/dl

Traceability:

GC-IDMS

Stability:

Without opening at 2-8°C until the expiry date.
After first opening at least 4 weeks at 2-8°C.


Warnings:


For in vitro diagnostic use only.

Do not use after the expiry date!

The following symbols are used on labels:

 For in vitro diagnostic use

 Expiry date (last day of the month)

 Temperature limitation

 Batch Code

 Code