

Gyártási szám: 7665 #####

Lejáratidő: ÉÉÉÉ-HH

Felhasználás:

Karbamid reagens kalibrálására szolgáló stabil folyékony karbamid standard.

Felhasználásra kész.

A standard koncentrációja LOT specifikus.

Koncentráció:

| Reagens | Érték | Egység |
|---|-------|--------|
| Karbamid stabil folyékony reagens (Ureáz/GLDH UV) | ## | mmol/l |



Visszavezethetőség:

SRM 909b

Stabilitás:Felbontás nélkül 2-8°C-on tárolva a lejáratidőig.
Felbontás után 2-8°C-on tárolva minimum 4 hétig.**Figyelmeztetések:**

Csak in vitro diagnosztikai célokra használható.

Ne használja a standardet a lejáratidő eltelte után!

A címkén a következő szimbólumok találhatóak: In vitro diagnosztikai használatra Lejáratidő (a hónap utolsó napja) Tárolási hőmérséklet Gyártási szám Kódszám

Lot No: 7665 #####

Expiry date: YYYY-MM

Preparation:

Urea standard for the calibration of urea reagent.

Ready to use.

Concentration is LOT specific.

Concentration:

| Reagent | Value | Unit |
|---|-------|-------|
| Urea stable liquid reagent (Urease/GLDH UV) | ## | mg/dl |



Traceability:

SRM 909b

Stability:Without opening at 2-8°C until the expiry date.
After first opening at least 4 weeks at 2-8°C.**Warnings:**

For in vitro diagnostic use only.

Do not use after the expiry date!

The following symbols are used on labels: For in vitro diagnostic use Expiry date (last day of the month) Temperature limitation Batch Code Code