

MAGNÉZIUM STABIL FOLYÉKONY
STANDARD



MAGNESIUM STABLE LIQUID STANDARD

Kódszám/Catalog No.: 52201
Kiszerezés/Size: 1x5 ml

Gyártási szám: 7185 #####

Lejárató idő: ÉÉÉÉ-HH

Felhasználás:

Magnézium reagens kalibrálására szolgáló stabil folyékony magnézium standard.
Felhasználásra kész.
A standard koncentrációja LOT specifikus.

Koncentráció:

| Reagens | Érték | Egység |
|-------------------------------------------------|-------|--------|
| Magnézium stabil folyékony reagens (xylidilkék) | ### | mmol/l |

Visszavezethetőség:

AAS


Stabilitás:


Felbontás nélkül 2-8°C-on tárolva a lejárató időig.
Felbontás után 2-8°C-on tárolva minimum 4 hétig.


Figyelmeztetések:

Csak in vitro diagnosztikai célokra használható.
Ne használja a standardet a lejárató dátum eltelte után!

A címkén a következő szimbólumok találhatóak:

 In vitro diagnosztikai használatra

 Lejárató idő (a hónap utolsó napja)

 Tárolási hőmérséklet

 Gyártási szám

 Kódszám

Lot No: 7185 #####

Expiry date: YYYY-MM

Preparation:

Magnesium standard for the calibration of magnesium reagent.
Ready to use.
Concentration is LOT specific.

Concentration:

| Reagent | Value | Unit |
|------------------------------------------------|-------|-------|
| Magnesium stable liquid reagent (xylydyl blue) | ### | mg/dl |

Traceability:

AAS


Stability:


Without opening at 2-8°C until the expiry date.
After first opening at least 4 weeks at 2-8°C.

Warnings:

For in vitro diagnostic use only.
Do not use after the expiry date!

The following symbols are used on labels:

 For in vitro diagnostic use

 Expiry date (last day of the month)

 Temperature limitation

 Batch Code

 Code